

Supplementary Material

Supplementary Figures

Figure 1 – CAD Model of 3D Clinostat v2 with sample holders

Figure 2 – Quantitating the Distribution of the Acceleration Vector

Figure 3 – Measurement of Incubator Temperature during Clinostat Operation

Figure 4 - Path and Distribution of Acceleration Vector for the 3D clinostat and RPM 2.0 when operated at inner frame velocity 1.5 rpm and outer frame velocity 3.875 rpm.

Figure 5 – Image of bacteria in flaskette after 4 days at normal gravity or simulated microgravity

Figure 6 – Venn diagram of comparisons of differential expression analyses

Figure 7 – Twenty-Four-Hour *M. marinum* Growth and Rifampicin Survival with data shown in OD₆₀₀, CFU/mL and % Rifampicin Survival

Figure 8 – Four-Day *M. marinum* Growth and Rifampicin Survival with data shown in OD₆₀₀, CFU/mL and % Rifampicin Survival

Figure 9 – Four-day *M. marinum* growth study

Supplementary Table

Table 1 – Number of genes with altered expression due to simulated microgravity

Supplementary Programs

Program 1 – Clinostat Control (Arduino C)

Program 2 – Arduino UNO Accelerometer (Arduino C)

Program 3 – Accelerometer Web Socket (Java Processing)

Program 4 – Computer Model (Python)

Program 5 – Data Processing and Visualization (Python)

Program 6 – Accelerometer Vector Path Visualization (Python)

Program 7 – Heatmap Generation (Python)

Program 8 – Temperature Sensor (Python)

Supplementary Figure 1




Figure 1 - CAD Model of 3D Clinostat v2 with sample holders

Fusion 360 was used to generate the designs for the 3D clinostat v2 (A), and the sample holders (B-D) that slide into the inner frame. The sample holder shown on the clinostat (A) can carry a cell culture flask with an area of 25 cm² (T25) and a flaskette on each side. It is necessary to balance the weight on each side of the holder. The HARV holder can carry two HARVs (B) and air flow to the membrane on the lower side of the HARV is not restricted due to the design of the holder. The T25 holder (C) can carry a T25 cell culture flask on each side of the holder, and the flaskette holder (D) can carry two flaskettes on each side. The samples are attached with zip ties using the holes in the holder.

Supplementary Figure 2




Figure 2 – Quantitating the Distribution of the Acceleration Vector

A Python program graphed the x, y, and z components of the acceleration vector and a program was developed to score how distributed the gravity path was around the unit sphere. The program evenly distributed 1000 points (yellow dots) around the unit sphere using a Fibonacci Lattice. For each acceleration vector position on the sphere (provided by the x, y and z), the program found the closest three points (yellow dots), effectively finding what triangular region the vector data point was located in (A). If the vector data point was the first to be located within that region, the distribution score was increased by one, otherwise, the score remained the same.

Supplementary Figure 3




Figure 3 - Incubator Temperature Change Due to Clinostat Operation

The 3D clinostat v1 was placed in a 37°C cell culture incubator. The temperature probe attached to a Raspberry Pi zero and a Pi Sugar battery (A) was placed in the incubator, and the incubator allowed to reach operating temperature. The temperature was monitored for 48 hours with either the 3D clinostat inactive (Control, B) or with the 3D clinostat operating with inner and outer rings at 0.625 rpm and 1.25 rpm, respectively. The average temperature was calculated from the data collected. The experiment was repeated 3 times, and the average temperature and the standard deviation calculated (B).

Supplementary Figure 4




Figure 4 – Path and Distribution of Acceleration Vector

An Arduino UNO WIFI Rev 2, which has an onboard accelerometer, was attached to the 3D clinostat and the RPM 2.0. The acceleration vector path was plotted using data obtained from the Arduino accelerometer with the machines operating at an inner frame and outer frame velocity combination of 1.5 rpm and 3.875 rpm, respectively (A and B). The data output by the algorithm from the RPM 2.0 was also plotted (C). The distribution score was calculated for these frame velocities and is shown in D. Data collected from 0-3000 seconds was used to calculate the distribution score. Path images (A-C) were from data collected between 2000-3000 seconds so the pattern could be seen. The conditions were performed three times and the average and standard deviation shown.

Supplementary Figure 5




Figure 5 – *M. marinum* form a mass of cells when exposed to simulated microgravity using the 3D clinostat v2 or RPM 2.0

Flaskettes were filled with bacteria in biofilm-promoting medium and attached to the RPM 2.0 or 3D clinostat. The RPM 2.0 was operated using the partial gravity mode and 0g file, and the 3D clinostat was operated with an inner frame velocity of 1.5 rpm and outer frame velocity of 3.875 rpm. Control flasks for normal gravity (NG) were placed on a shelf in the same incubator that housed the 3D clinostat and RPM 2.0. After four days at 31°C, the flaskettes were removed and the images captured. The bacteria were in a mass when removed from the 3D clinostat and the RPM 2.0, but were dispersed aggregates in the NG flask.

Supplementary Figure 6




Figure 6 – Venn diagram showing the overlap of gene expression changes induced by growth of *M. marinum* when exposed to simulated microgravity using the RPM 2.0 and 3D clinostat

Differential expression analysis for 3071 genes was performed on RNA-Seq data from RNA isolated from *M. marinum* grown in biofilm-promoting medium for 4 days under normal gravity (NG) or simulated microgravity on the RPM 2.0 or the 3D clinostat (I: 1.5 rpm O: 3.875 rpm). Comparisons of the differential analyses for RPM 2.0 and NG, and the 3D Clinostat and NG demonstrated an overlap of 426 genes.

Supplementary Figure 7




Figure 7 –Twenty-Four-Hour *M. marinum* Growth and Rifampicin Survival

M. marinum were grown in 10 mL HARVs in medium to promote planktonic growth with or without 0.8 µg/mL rifampicin for 24 hrs at 31°C. The RPM 2.0 (RPM) was operated in partial gravity mode using the 0 g file, and the 2D clinostat from Synthecon was used for the normal gravity (NG) and LSMMG conditions at 25 rpm. For the 3D clinostat, I = inner frame velocity in rpm and O = outer frame velocity in rpm. Work was performed at Kennedy Space Center in A, C and E and at LSUHS for B, D and F. The OD₆₀₀ (A and B) and the CFU/mL (C and D) were used to compare growth without rifampicin. The percentage of *M. marinum* surviving the 24-hr rifampicin treatment is shown in D and E. Experiments were performed three times and the average and standard deviation is shown. * p < 0.05 compared to NG, # p < 0.05 compared to RPM, \$ p < 0.05 compared to LSMMG

Supplementary Figure 8




Figure 8 –Four-Day *M. marinum* Growth and Rifampicin Survival

M. marinum were grown in 10 mL HARVs in medium to promote planktonic growth with or without 0.05 µg/mL rifampicin for 4 days at 31°C. The 2D clinostat from Synthecon operated at 25 rpm for the NG and LSMMG conditions. For the 3D clinostat, I = inner frame velocity in rpm and O = outer frame velocity in rpm. All work was performed at LSUHS. Experiments in A, C and E used the 4 unit Synthecon 2D clinostat for the NG and LSMMG cultures. Due to the size of the units, the NG and LSMMG cultures were grown in a different incubator to the 3D clinostat culture. Experiments in B, D and F used two single unit Synthecon 2D clinostats for LSMMG and the cultures were grown in the same incubator as the 3D clinostat culture. The OD₆₀₀ (A and B) and CFU/ mL (C and D) are shown for cultures without rifampicin. The percentage of *M. marinum* surviving the 4-day rifampicin treatment is shown in E and F. Experiments were performed three times and the average and standard deviation is shown. * p < 0.05 compared to NG, \$ p < 0.05 compared to LSMMG, # p < 0.05 compared to I:2 rpm O: 1.125 rpm

Supplementary Figure 9




Figure 9 –Four day *M. marinum* growth

M. marinum were grown in 10 mL HARVs in medium to promote planktonic growth for 4 days at 31°C. The 2D clinostat from Synthecon was operated at 25 rpm for the normal gravity (NG) and LSMMG conditions. For the 3D clinostat, I = inner frame velocity in rpm and O = outer frame velocity in rpm. Experiments were performed three times and the average and standard deviation is shown for the OD₆₀₀ (A) and the CFU/ mL (B). The p values are in comparison to the 2D LSMMG culture.

Table 1 - Number of genes with altered transcript levels due to simulated microgravity

The rRNA genes were removed from the RNA-seq data before performing the analyses below. Genes with low expression were removed if the minimum count was less than 10 (**A** and **B**) or 5 (**C** and **D**) in approximately 70% of the samples, and the minimum total count of all samples had to equal at least 15. In **E**, low expression genes were not filtered from the analysis.

A: The outlier 3D clinostat sample was included in the analysis. Minimum count less than 10.

Type of Gene expression change	Type of Comparison		
	3D clinostat v. RPM 2.0	RPM 2.0 v. NG	3D clinostat v. NG
Decreased expression	0	350	292
No change	3071	2474	2597
Increased expression	0	247	182

B: The outlier 3D clinostat sample was excluded from the analysis. Minimum count less than 10.

Type of Gene expression change	Type of Comparison		
	3D clinostat v. RPM 2.0	RPM 2.0 v. NG	3D clinostat v. NG
Decreased expression	0	393	319
No change	3071	2405	2530
Increased expression	0	273	222

C: The outlier 3D clinostat sample was included in the analysis. Minimum count less than 5.

Type of Gene expression change	Type of Comparison		
	3D clinostat v. RPM 2.0	RPM 2.0 v. NG	3D clinostat v. NG
Decreased expression	0	383	320
No change	3925	3281	3414
Increased expression	0	261	191

D: The outlier 3D clinostat sample was excluded from the analysis. Minimum count less than 5.

Type of Gene expression change	Type of Comparison		
	3D clinostat v. RPM 2.0	RPM 2.0 v. NG	3D clinostat v. NG
Decreased expression	0	428	332
No change	3925	3216	3385
Increased expression	0	281	208

E: The outlier 3D clinostat sample was excluded from the analysis. No filtering for low gene expression.

Type of Gene expression change	Type of Comparison		
	3D clinostat v. RPM 2.0	RPM 2.0 v. NG	3D clinostat v. NG
Decreased expression	0	480	374
No change	5557	4839	5016
Increased expression	0	238	167

Supplementary Program 1 – Clinostat Control (Arduino C)

```
1 #include <ax12.h>
2 #include <BioloidController.h>
3
4 void setup()
5 {
6     //Set motor to wheel mode
7     ax12SetRegister2(3, AX_CCW_ANGLE_LIMIT_L, 0); //inner
8     ax12SetRegister2(2, AX_CCW_ANGLE_LIMIT_L, 0); //outer
9     delay(1000);
10 }
11
12 void loop()
13 {
14     ax12SetRegister2(2, AX_GOAL_SPEED_L, 50); //inner
15     ax12SetRegister2(3, AX_GOAL_SPEED_L, 4); //outer
16     Serial.println(inner);
17     delay(250);
18 }
19
```

Supplementary Program 2 – Arduino UNO Accelerometer (Arduino C)

```
1 #include <ArduinoHttpClient.h>
2 #include <WiFiNINA.h>
3 #include <Arduino_LSM6DS3.h>
4 #include "arduino_secrets.h"
5
6 //Set Wifi credentials in arduino_secrets.h
7 //Wifi Settings
8 char ssid[] = SECRET_SSID;
9 char pass[] = SECRET_PASS;
10
11 char serverAddress[] = "192.168.1.1"; //SET RECEIVING COMPUTER IP ADDRESS
12 int port = 8080;
13
14 WiFiClient wifi;
15 WebSocketClient client = WebSocketClient(wifi, serverAddress, port);
16 int status = WL_IDLE_STATUS;
17
18 int failCount = 0;
19
20 unsigned long t;
21
22
23 void setup() {
24   Serial.begin(9600);
25   wifi_setup();
26
27   //Setup for Accelerometer
28   if (!IMU.begin()){
29     Serial.println("Failed to initialize IMU!");
30     while (1);
31   }
32
33   Serial.print("Accelerometer sample rate = ");
34   Serial.print(IMU.accelerationSampleRate());
35   Serial.println(" Hz");
36   Serial.println();
37   Serial.println("Acceleration in G's");
38   Serial.println("X\tY\tZ");
39 }
40
41 void wifi_setup() {
42   //Setup for WIFI
43   while (status != WL_CONNECTED) {
44     Serial.println(status);
45     Serial.print("Attempting to connect to Network named: ");
46     Serial.println(ssid);           // print the network name (SSID);
47
48     // Connect to WPA/WPA2 network:
49     status = WiFi.begin(ssid, pass);
50   }
51
52   // print the SSID of the network you're attached to:
53   Serial.print("SSID: ");
54   Serial.println(WiFi.SSID());
55
56   // print your WiFi shield's IP address:
57   IPAddress ip = WiFi.localIP();
58   Serial.print("IP Address: ");
59   Serial.println(ip);
```

```

60 }
61
62 void loop() {
63     float x,y,z;
64     int messageSize;
65
66     Serial.println("Starting WebSocket client");
67     client.begin();
68     Serial.println("Websocket Connected");
69
70     while (client.connected()) {
71         t = millis();
72         if (IMU.accelerationAvailable()) {
73             IMU.readAcceleration(x,y,z);
74
75             client.beginMessage(TYPE_TEXT);
76             client.print(x);
77             client.print('\t');
78             client.print(y);
79             client.print('\t');
80             client.print(z);
81             client.print('\t');
82             client.println(t);
83             client.endMessage();
84
85             Serial.print(x);
86             Serial.print('\t');
87             Serial.print(y);
88             Serial.print('\t');
89             Serial.println(z);
90             messageSize = client.parseMessage();
91
92             if (messageSize == 0) {
93                 failCount++;
94                 Serial.print("Fail: ");
95                 Serial.println(failCount);
96             }
97             else {
98                 failCount = 0;
99             }
100
101             if (failCount > 10) {
102                 break;
103             }
104         }
105     }
106     delay(1000);
107 }
108 client.stop();
109 wifi_setup();
110 client.println("disconnected\n");
111 Serial.println("disconnected\n");
112 }
113

```

Supplementary Program 3 - Accelerometer Web Socket (Java Processing)

```
1 import websockets.*;
2
3 WebSocketServer ws;
4 int now;
5 int count = 0;
6 int old_count = 0;
7 int failCount = 0;
8 float x, y;
9 String old_msg;
10 PrintWriter output;
11
12 String path = "C:/OUTPATH/"; //SET OUTPUT PATH HERE
13 String name = path + "FILENAME.txt"; //SET OUTPUT FILE NAME HERE
14
15 void setup() {
16     size(200, 200);
17     ws = new WebSocketServer(this, 8080, "");
18
19     output = createWriter(name);
20     now=millis();
21     x=0;
22     y=0;
23 }
24
25 void draw() {
26     //Debug display
27     ellipse(x, y, 10, 10);
28
29     if (millis()>now+1000) {
30         if(old_count == count) {
31             print("Fail ");
32             print(failCount);
33             print("/");
34             println(count);
35             output.print("fail: ");
36             output.println(millis());
37             output.flush();
38             background(random(255), random(255), random(255));
39             failCount++;
40         }
41         else {
42             ws.sendMessage("Success");
43             old_count = count;
44         }
45
46         now = millis();
47     }
48 }
49
50 void websocketServerEvent(String msg) {
51     //Triggered when data received
52     //Outputs data to file that can be parsed by python program
53     count++;
54     println(msg);
55     output.println(msg);
56     output.flush();
57     x=random(width);
58     y=random(height);
59 }
```

Supplementary Program 4 – Computer Model (Python)

```
1 import math
2
3 class Sim:
4     #Returns x-component of the G vector at a specific time for a specific inner and
    #outer frame angular velocities
5     def gVectorX(self, timeInSeconds, localInnerInRadSec, localOuterInRadSec):
6         ret =
7         math.sin(localOuterInRadSec*timeInSeconds)*math.cos(localInnerInRadSec*timeInSeconds)
8         return ret
9     #Returns z-component of the G vector at a specific time for a specific outer
    #frame angular velocity
10    def gVectorY(self, timeInSeconds, localOuterInRadSec):
11        ret = math.cos(localOuterInRadSec*timeInSeconds)
12        return ret
13    #Returns y-component of the G vector at a specific time for a specific inner and
    #outer frame angular velocities
14    def gVectorZ(self, timeInSeconds, localInnerInRadSec, localOuterInRadSec):
15        ret =
16        math.sin(localOuterInRadSec*timeInSeconds)*math.sin(localInnerInRadSec*timeInSeconds)
17        return ret
18
19    #Converts and returns frame velocity from RPM to radians per second
20    def RPMtoRadSec(self, RPM):
21        ret = RPM * (math.pi / 30)
22        return ret
23
24    #Returns a Python dictionary containing value of g vector of each axis at
    #specific times
25    def gVectorData(self, startTimeInSeconds, endTimeInSeconds, innerRPM, outerRPM):
26        innerInRadSec = self.RPMtoRadSec(innerRPM)
27        outerInRadSec = self.RPMtoRadSec(outerRPM)
28        timeArray, xArray, yArray, zArray = [], [], [], []
29        for t in range(startTimeInSeconds, endTimeInSeconds + 1):
30            timeArray.append(t)
31            xArray.append(self.gVectorX(t, innerInRadSec, outerInRadSec))
32            yArray.append(self.gVectorY(t, outerInRadSec))
33            zArray.append(self.gVectorZ(t, innerInRadSec, outerInRadSec))
34        data = timeArray, xArray, yArray, zArray
35        return data
```


Supplementary Program 5 - Data Processing and Visualization (Python)

```
1 import os
2 import sys
3 import numpy as np
4 from pathVisualization import PathVisualization
5 from simulation import Sim
6
7 class DataProcessor:
8     def __init__(self):
9         #Uncomment machine type based on data source
10         self.machineType = "sim"
11         #self.machineType = "clinostat"
12         #self.machineType = "rpm-integ"
13         #self.machineType = "rpm-ard"
14
15         #Set range for average magnitude (seconds)
16         self.minSeg = 2000
17         self.maxSeg = 3000
18
19         #Set data output path
20         outParent = "C:\\\\OUTPATH\\" #SET OUTPUT FILE PATH
21
22         if self.machineType != "sim":
23             self.inPath = "C:\\\\INPATH\\" #SET INPUT FILE PATH
24             self.inName = "exampleFile.txt" #SET FILE NAME
25
26             try:
27                 testOpen = open(self.inPath + self.inName, "r")
28                 testOpen.close()
29             except FileNotFoundError:
30                 print("\nERROR: Input file not found: " + self.inPath + self.inName +
31 "\n")
32                 sys.exit()
33             if self.machineType != "sim":
34                 self.outPath = outParent + self.inName[:len(self.inName)-4] + "\\\"
35             else:
36                 self.outPath = outParent + self.inName + "\\\"
37             print(self.outPath)
38         else: #Applies only for Computer Model
39             self.innerV = 3.875 #SET INNER VELOCITY
40             self.outerV = 1.5 #SET OUTER VELOCITY
41             self.inName = str(self.innerV) + "-" + str(self.outerV)
42             self.outPath = outParent + self.inName + "\\\"
43
44             if os.path.isdir(self.outPath) == False:
45                 os.makedirs(self.outPath)
46
47     def outputDataFile(self):
48         """Outputs processed data file"""
49         self.__getVectors()
50         xTimeAvg, yTimeAvg, zTimeAvg = self.__getTimeAvg()
51         magList = self.__getMagnitude(xTimeAvg, yTimeAvg, zTimeAvg)
52         avgMagSeg = self.__getMagSeg(magList)
53         disScore = self.__getDisScore()
54
55         fileOut = open(self.outPath + 'processedAccelData.txt', 'w+')
56
57         if self.machineType != 'sim':
58             fileOut.write(self.inName)
59         else:
```

```

59         fileOut.write('Inner: ' + str(self.innerV) + '\t' + 'Outer: ' +
60         str(self.outerV) + '\n')
61
62         fileOut.write('Segment: ' + str(self.minSeg) + ":" + str(self.maxSeg) + '\t'
63         + str(avgMagSeg) + '\t' + '\t\t\t' + 'Dispersion: ' + '\t' + str(disScore) + '\n')
64         fileOut.write('Time' + '\t' + 'X Vector' + '\t' + 'Y Vector' + '\t' + 'Z
65         Vector' + '\t' + 'X Time Avg' + '\t' + 'Y Time Avg' + '\t' + 'Z Time Avg' + '\t' +
66         'Magnitude' + '\n')
67
68         for i in range(len(self.x)):
69             fileOut.write(str(self.time[i]) + '\t' + str(self.x[i]) + '\t' +
70             str(self.y[i]) + '\t' + str(self.z[i]) + '\t' + str(xTimeAvg[i]) + '\t' +
71             str(yTimeAvg[i]) + '\t' + str(zTimeAvg[i]) + '\t' + str(magList[i]) + '\n')
72
73         fileOut.close()
74
75         print("Average Magnitude: " + str(avgMagSeg))
76         print("Distribution: " + str(disScore))
77
78     def __getVectors(self):
79         """Directs where to get data from"""
80         if self.machineType == 'sim':
81             self.time, self.x, self.y, self.z = self.__getSimAccelData()
82         elif self.machineType == 'clinostat' or self.machineType == 'rpm-ard':
83             self.time, self.x, self.y, self.z = self.__getArdAccelData()
84         else:
85             self.time, self.x, self.y, self.z = self.__getRPMAccelData()
86
87     def __getArdAccelData(self):
88         """Reads data from file formatted by the Arduino UNO Accelerometer"""
89         dataIn = open(self.inPath + self.inName, 'r')
90
91         x, y, z, time = [], [], [], []
92         for line in dataIn.readlines():
93             lineSplit = line.split('\t')
94             if len(lineSplit) == 4:
95                 x.append(float(lineSplit[0]))
96                 y.append(float(lineSplit[1]))
97                 z.append(float(lineSplit[2]))
98                 unTime = float(lineSplit[3]) / 1000
99                 time.append(unTime)
100
101         dataIn.close()
102         return time, x, y, z
103
104     def __getRPMAccelData(self):
105         """Reads data from file formatted by the RPM 2.0"""
106         dataIn = open(self.inPath + self.inName, 'r')
107
108         x, y, z, time = [], [], [], []
109         for line in dataIn.readlines():
110             lineFirstSplit = line.split('\t')
111             splitVectors = lineFirstSplit[1]
112             vectorList = splitVectors.split(' ')
113
114             x.append(float(vectorList[0]))
115             z.append(float(vectorList[1]))
116             y.append(float(vectorList[2]))
117
118         for i in range(len(x)):

```

```

113         time.append(i)
114
115     dataIn.close()
116     return time, x, y, z
117
118     def __getSimAccelData(self):
119         """Gets data from computer model"""
120         simInnerV = float(self.innerV)
121         simOuterV = float(self.outerV)
122         vectorSim = Sim()
123         time, x, y, z = vectorSim.gVectorData(0, self.maxSeg, simInnerV, simOuterV)
124         return time, x, y, z
125
126     def __getTimeAvg(self):
127         """Calculates Time Average"""
128         xTimeAvg = []
129         xTempList = []
130         for xIter in self.x:
131             xTempList.append(xIter)
132             xTimeAvg.append(np.mean(xTempList))
133
134         yTimeAvg = []
135         yTempList = []
136         for yIter in self.y:
137             yTempList.append(yIter)
138             yTimeAvg.append(np.mean(yTempList))
139
140         zTimeAvg = []
141         zTempList = []
142         for zIter in self.z:
143             zTempList.append(zIter)
144             zTimeAvg.append(np.mean(zTempList))
145
146         return xTimeAvg, yTimeAvg, zTimeAvg
147
148     def __getMagnitude(self, xTimeAvg, yTimeAvg, zTimeAvg):
149         """Calculates magnitude of time-averaged acceleration vector"""
150         magList = []
151
152         for i in range(len(self.x)):
153             xIter = xTimeAvg[i]
154             yIter = yTimeAvg[i]
155             zIter = zTimeAvg[i]
156
157             mag = (xIter ** 2 + yIter ** 2 + zIter ** 2) ** 0.5
158             magList.append(mag)
159
160         return magList
161
162     def __getMagSeg(self, magList):
163         """Calculates average magnitude between segments"""
164         magSegList = magList[self.minSeg:self.maxSeg]
165         if len(magList) < self.minSeg:
166             print("\nERROR: Segment begins after data ends - " + str(len(magList)) +
167 " sec\n")
168             sys.exit()
169         elif len(magSegList) < (self.maxSeg - self.minSeg):
170             print("\nWARNING: Not enough data for segment - " + str(len(magList)) + "
171 sec\n")
172         return np.mean(magList[self.minSeg:self.maxSeg])

```



```

171
172 def __getDisScore(self):
173     xSeg = self.x[:self.maxSeg]
174     ySeg = self.y[:self.maxSeg]
175     zSeg = self.z[:self.maxSeg]
176     path = PathVisualization(self.inName, xSeg, ySeg, zSeg)
177     disScore = path.getDistribution()
178
179     return disScore
180
181 def copyRawData(self):
182     if self.machineType == "sim":
183         print("\nWARNING: Machine Type = Sim - No files were copied\n")
184         return
185     dataIn = open(self.inPath + self.inName, 'r')
186     copyData = open(self.outPath + self.inName[:len(self.inName)-4] + '-Copy' +
self.inName[len(self.inName)-4:], 'w+')
187
188     copyData.write(self.inPath + self.inName + '\n\n')
189
190     for line in dataIn.readlines():
191         copyData.write(line)
192
193     copyData.close()
194     dataIn.close()
195
196 def createGraphs(self):
197     self.__getVectors()
198     xTimeAvg, yTimeAvg, zTimeAvg = self.__getTimeAvg()
199     magnitude = self.__getMagnitude(xTimeAvg, yTimeAvg, zTimeAvg)
200
201     if self.machineType != 'sim':
202         graphPath = PathVisualization(self.inName, self.x, self.y, self.z,
saveFile=self.outPath + self.inName[:len(self.inName)-4])
203     else:
204         graphPath = PathVisualization(self.inName, self.x, self.y, self.z,
saveFile=self.outPath + self.inName)
205
206     graphPath.createPathShadowFig(mode='save', legend=False)
207
208     graphPath.createVectorFig(self.time, mode='save')
209     graphPath.createTimeAvgFig(xTimeAvg, yTimeAvg, zTimeAvg, self.time,
mode='save')
210     graphPath.createMagFig(magnitude, self.time, mode='save')
211
212 def createSegGraphs(self):
213     self.__getVectors()
214     xSeg = self.x[:self.maxSeg]
215     ySeg = self.y[:self.maxSeg]
216     zSeg = self.z[:self.maxSeg]
217
218     if self.machineType != 'sim':
219         comparePaths = PathVisualization(self.inName, xSeg, ySeg, zSeg,
saveFile=self.outPath + self.inName[:len(self.inName)-4] + '-Segment')
220     else:
221         comparePaths = PathVisualization(self.inName, xSeg, ySeg, zSeg,
saveFile=self.outPath + self.inName + '-Segment')
222     comparePaths.createPathFig(mode='save')
223
224 if __name__ == "__main__":

```

```
225 process = DataProcessor()  
226 process.outputDataFile()  
227 process.createGraphs()  
228 process.createSegGraphs()  
229 process.copyRawData()  
230  
231
```

Supplementary Program 6 - Accelerometer Vector Path Visualization (Python)

```
1 import numpy as np
2 import matplotlib.pyplot as plt
3 import matplotlib.cm as cm
4 from mpl_toolkits.mplot3d import Axes3D
5
6 class PathVisualization:
7     def __init__(self, ID, x, y, z, saveFile=''):
8         self.ID = ID
9
10        self.x = x
11        self.y = y
12        self.z = z
13
14        self.pathCoords = list(zip(self.x, self.y, self.z))
15        self.num_points = 1000
16
17        self.saveFile = saveFile
18
19    def __createSphere(self):
20        """creates guide sphere with fibonnaci sequence on same plot as path"""
21        golden_r = (np.sqrt(5.0) + 1.0) / 2.0      #golden ratio = 1.61803...
22        golden_a = (2.0 - golden_r) * (2.0 * np.pi)  #golden angle + 2.39996...
23
24        Xs, Ys, Zs = [], [], []
25
26        for i in range(self.num_points):
27            ys = 1 - (i / float(self.num_points - 1)) * 2
28            radius = np.sqrt(1 - ys * ys)
29
30            theta = golden_a * i
31
32            xs = np.cos(theta) * radius
33            zs = np.sin(theta) * radius
34
35            Xs.append(xs)
36            Ys.append(ys)
37            Zs.append(zs)
38
39        return(Xs, Ys, Zs)
40
41    def __splitSphere(self, sphereCoords):
42        octants = {'posI':[], 'posII':[], 'posIII':[], 'posIV':[], 'negI':[],
43        'negII':[], 'negIII':[], 'negIV':[]}
44
45        for row in sphereCoords:
46            #Coordinates: x = row[0], y = row[1], z = row[2]
47            if (row[2] > 0):
48                if (row[1] > 0):
49                    if (row[0] > 0):
50                        octants['posI'].append(row)
51                    else:
52                        octants['posII'].append(row)
53                elif (row[0] > 0):
54                    octants['posIV'].append(row)
55                else:
56                    octants['posIII'].append(row)
57            else:
58                if (row[1] > 0):
59                    if (row[0] > 0):
```

```

59         octants['negI'].append(row)
60     else:
61         octants['negII'].append(row)
62 elif (row[0] > 0):
63     octants['negIV'].append(row)
64 else:
65     octants['negIII'].append(row)
66
67     return(octants)
68
69 def __getPathOctant(self, pathRow):
70     """Checks what quadrant the tuple of coordinates is in"""
71     if (pathRow[2] > 0):
72         if (pathRow[1] > 0):
73             if (pathRow[0] > 0):
74                 return 'posI'
75             else:
76                 return 'posII'
77         elif (pathRow[0] > 0):
78             return 'posIV'
79         else:
80             return 'posIII'
81     else:
82         if (pathRow[1] > 0):
83             if (pathRow[0] > 0):
84                 return 'negI'
85             else:
86                 return 'negII'
87         elif (pathRow[0] > 0):
88             return 'negIV'
89         else:
90             return 'negIII'
91
92 def __getDistanceBetween(self, pathTupleCoords, sphereTupleCoords):
93     """Calculates distance between two points (path point and sphere point)"""
94     pathX, pathY, pathZ = pathTupleCoords
95     sphereX, sphereY, sphereZ = sphereTupleCoords
96
97     diffX = pathX - sphereX
98     diffY = pathY - sphereY
99     diffZ = pathZ - sphereZ
100
101     sumSquares = diffX ** 2 + diffY ** 2 + diffZ ** 2
102     dist = np.sqrt(sumSquares)
103
104     return dist
105
106 def __getDistributionNum(self, sphereCoords):
107     """Creates a dictionary of the three sphere coords that are closest to each
108     path coord.
109     The length of this dictionary will be the distribution score assigned to the
110     combination.
111     A higher distribution score means the path of the combination will better
112     cover the entire sphere,
113     thus the average orientation will be closer to zero and create a better
114     simulation."""
115     octants = self.__splitSphere(sphereCoords) #splits spheres into octants and
116     returns dictionary with each octant separated

```



```

113     #Every time the path passes between a new segment, the three vertices are
114     added as keys.
115     #The value is a list of the indexes of the pathCoords to denote the how long
116     it takes to repeat the entire path
117     pathMap = {}
118     repeatTime = -1
119
120     for pathRow in self.pathCoords:
121         pathOctant = self.__getPathOctant(pathRow) #determines what octant this
122         set of path coordinates is in
123         sphereCoordsSplit = octants[pathOctant] #retrieves sphere coordinates in
124         the octant that this set of path coordinates is in
125         distDict = {}
126         repeatTime += 1
127
128         for sphereRow in sphereCoordsSplit:
129             dist = self.__getDistanceBetween(pathRow, sphereRow) #gets distance
130             between two points
131             distDict[sphereRow] = dist #Updates dictionary (key:sphere
132             coordinates, value: distance between path coordinate and that sphere coordinate)
133
134             rankedDist = sorted(distDict.items(), key=lambda x:x[1]) #takes
135             dictionary of distances, sorts them, and makes it into a list of tuples
136             segmentVertices = (rankedDist[0][0], rankedDist[1][0], rankedDist[2][0])
137             #the three closest sphere coordinates are the vertices of the segment that the path
138             coordinate is in
139
140             pathMap[segmentVertices] = pathMap.get(segmentVertices, []) +
141             [repeatTime] # If the path has not passed through this segment, it creates a key and
142             assigns the value as the first time the path passed through the segment
143
144         return(len(pathMap))
145
146     def getDistribution(self):
147         #get sphere coords
148         Xsphere, Ysphere, Zsphere = self.__createSphere()
149         sphereCoords = list(zip(Xsphere, Ysphere, Zsphere))
150         score = self.__getDistributionNum(sphereCoords)
151         return score
152
153     """
154     Distribution Figure Types:
155     Distribution      --> Shadow
156     createPathFig     --> Path
157     createSphereFig   --> Points
158     createPathShadowFig --> Path, Shadow
159     createCombinedFig --> Path, Shadow, Points
160     """
161
162     def createPathShadowFig(self, mode='save', separated = False, title = True,
163     legend = True):
164         disScore = self.getDistribution()
165         legTitle = 'Distribution: ' + str(disScore)
166         if separated:
167             fig = plt.figure(figsize=plt.figaspect(0.4)) #Separates Shadow and Path
168             plt.xlim(-1.0, 1.0)
169             plt.ylim(-1.0, 1.0)
170
171             if title:
172                 fig.suptitle(self.ID + ' -- ' + 'Path and Shadow')
173
174

```



```

161         ax = fig.add_subplot(1, 2, 1, projection='3d')
162         ax.plot(self.x, self.y, self.z, label = "Acceleration Vector Path",
linewidth = 1) #Plots Path
163         ax.legend()
164
165         ax = fig.add_subplot(1, 2, 2, projection='3d')
166         ax.plot(self.x, self.y, zs=-1, zdir='z', label='Projection in (x,y)')
#Plots Shadow on x,y plane
167         ax.plot(self.x, self.z, zs=1, zdir='y', label='Projection in (x,z)')
#Plots Shadow on x,z plane
168         ax.plot(self.y, self.z, zs=-1, zdir='x', label='Projection in (y,z)')
#Plots Shadow on y,z plane
169
170
171         if legend:
172             ax.legend(title=legTitle, loc='lower left')
173
174         else:
175             fig = plt.figure(figsize=plt.figaspect(0.85))
176             plt.xlim(-1.0, 1.0)
177             plt.ylim(-1.0, 1.0)
178
179             if title:
180                 fig.suptitle(self.ID + ' -- ' + 'Path and Shadow')
181
182         ax = fig.add_subplot(1, 1, 1, projection='3d') #Combines Shadow and Path
(default)
183         ax.plot(self.x, self.y, self.z, label = "Acceleration Vector Path",
linewidth = 1) #Plots Path
184         ax.plot(self.x, self.y, zs=-1, zdir='z', alpha = 0.3, label='Projection
in (x,y)') #Plots Shadow on x,y plane
185         ax.plot(self.x, self.z, zs=1, zdir='y', alpha = 0.3, label='Projection in
(x,z)') #Plots Shadow on x,z plane
186         ax.plot(self.y, self.z, zs=-1, zdir='x', alpha = 0.3, label='Projection
in (y,z)') #Plots Shadow on y,z plane
187
188         if legend:
189             ax.legend(title=legTitle, loc='lower left')
190
191         if mode == 'save':
192             plt.savefig(self.saveFile + '-path.png')
193         elif mode == 'show':
194             plt.show()
195         else:
196             plt.savefig(self.saveFile + '-path.png')
197             plt.show()
198
199         def createCombinedFig(self, mode='save', separated = False, title = True, legend
= True):
200             pointSize = np.ones(self.num_points) * 2
201             Xsphere, Ysphere, Zsphere = self.__createSphere()
202
203             if separated:
204                 fig = plt.figure(figsize=plt.figaspect(0.4)) #Separates Shadow and Path
205                 plt.xlim(-1.0, 1.0)
206                 plt.ylim(-1.0, 1.0)
207
208                 if title:
209                     fig.suptitle(self.ID + ' -- ' + 'Combined')
210

```

```

211         ax = fig.add_subplot(1, 2, 1, projection='3d')
212         ax.scatter(Xsphere, Ysphere, Zsphere, s = pointSize, color = 'orange')
        #Plots Sphere points
213         ax.plot(self.x, self.y, self.z, label = "Acceleration Vector Path",
linewidth = 1) #Plots Path
214
215         if legend:
216             ax.legend()
217
218         ax = fig.add_subplot(1, 2, 2, projection='3d')
219         ax.plot(self.x, self.y, zs=-1, zdir='z', label='Projection in (x,y)')
        #Plots Shadow on x,y plane
220         ax.plot(self.x, self.z, zs=1, zdir='y', label='Projection in (x,z)')
        #Plots Shadow on x,z plane
221         ax.plot(self.y, self.z, zs=-1, zdir='x', label='Projection in (y,z)')
        #Plots Shadow on y,z plane
222
223         if legend:
224             ax.legend()
225
226     else:
227         fig = plt.figure(figsize=plt.figaspect(0.85))
228         plt.xlim(-1.0, 1.0)
229         plt.ylim(-1.0, 1.0)
230
231         if title:
232             fig.suptitle(self.ID + ' -- ' + 'Combined')
233
234         ax = fig.add_subplot(1, 1, 1, projection='3d') #Combines Shadow and Path
        (default)
235         ax.scatter(Xsphere, Ysphere, Zsphere, s = pointSize, color = 'goldenrod')
        #Plots Sphere points
236         ax.plot(self.x, self.y, self.z, label = "Acceleration Vector Path",
linewidth = 1) #Plots Path
237         #ax.plot(self.x, self.y, zs=-1, zdir='z', alpha = 0.3, label='Projection
in (x,y)') #Plots Shadow on x,y plane
238         #ax.plot(self.x, self.z, zs=1, zdir='y', alpha = 0.3, label='Projection
in (x,z)') #Plots Shadow on x,z plane
239         #ax.plot(self.y, self.z, zs=-1, zdir='x', alpha = 0.3, label='Projection
in (y,z)') #Plots Shadow on y,z plane
240
241         if legend:
242             ax.legend()
243
244         if mode == 'save':
245             plt.savefig(self.saveFile + '-combinedPathSphere.png')
246         elif mode == 'show':
247             plt.show()
248         else:
249             plt.savefig(self.saveFile + '-combinedPathSphere.png')
250             plt.show()
251
252     def createSphereFig(self, mode='save', title = True):
253         fig = plt.figure(figsize=plt.figaspect(0.85))
254         plt.xlim(-1.0, 1.0)
255         plt.ylim(-1.0, 1.0)
256
257         if title:
258             fig.suptitle("Number of Points: " + str(self.num_points))
259

```

```

260     ax = fig.add_subplot(1, 1, 1, projection='3d')
261
262     pointSize = np.ones(self.num_points) * 2
263     Xsphere, Ysphere, Zsphere = self.__createSphere()
264     ax.scatter(Xsphere, Ysphere, Zsphere, s = pointSize, color = 'goldenrod')
265
266     if mode == 'save':
267         plt.savefig(self.saveFile + '-sphere.png')
268     elif mode == 'show':
269         plt.show()
270     else:
271         plt.savefig(self.saveFile + '-sphere.png')
272         plt.show()
273
274     def createShadowFig(self, mode='save', title = True, legend = True):
275         fig = plt.figure(figsize=plt.figaspect(0.85))
276         plt.xlim(-1.0, 1.0)
277         plt.ylim(-1.0, 1.0)
278
279         if title:
280             fig.suptitle(self.ID + ' -- ' + 'Shadow')
281
282         ax = fig.add_subplot(1, 1, 1, projection='3d')
283         ax.plot(self.x, self.y, zs=-1, zdir='z', label='Projection in (x,y)') #Plots
Shadow on x,y plane
284         ax.plot(self.x, self.z, zs=1, zdir='y', label='Projection in (x,z)') #Plots
Shadow on x,z plane
285         ax.plot(self.y, self.z, zs=-1, zdir='x', label='Projection in (y,z)') #Plots
Shadow on y,z plane
286
287         if legend:
288             ax.legend()
289
290         if mode == 'save':
291             plt.savefig(self.saveFile + '-shadowFig.png')
292         elif mode == 'show':
293             plt.show()
294         else:
295             plt.savefig(self.saveFile + '-shadowFig.png')
296             plt.show()
297
298     def createPathFig(self, mode='save', title = True):
299         fig = plt.figure(figsize=plt.figaspect(0.85))
300         plt.xlim(-1.0, 1.0)
301         plt.ylim(-1.0, 1.0)
302
303         if title:
304             fig.suptitle(self.ID + ' -- ' + 'Path')
305
306         ax = fig.add_subplot(1, 1, 1, projection='3d')
307         ax.plot(self.x, self.y, self.z, label = "Acceleration Vector Path", linewidth
= 1) #Plots Path
308
309         if mode == 'save':
310             plt.savefig(self.saveFile + '-pathFig.png')
311         elif mode == 'show':
312             plt.show()
313         else:
314             plt.savefig(self.saveFile + '-pathFig.png')
315             plt.show()

```



```

316
317
318 Acceleration Figure Types:
319 createVectorFig --> Figure with individual graphs of x, y, and z raw accel
values over time
320 createTimeAvgFig --> Figure with individual graphs of time averages of x, y, and
z over time
321 createMagFig --> Figure with magnitude over time
322
323
324 def formatTime(self, time):
325     startTime = time[0]
326     fTime = []
327     for t in time:
328         fTime.append(t - startTime)
329     return fTime
330
331 def createVectorFig(self, time, mode='save', title=True):
332     #fTime = self.formatTime(time)
333
334     #Plot X Vector
335     fig = plt.figure()
336     ax = fig.add_subplot(1, 1, 1)
337     ax.plot(time, self.x)
338
339     if title:
340         fig.suptitle(self.ID + ' -- ' + 'X Vector')
341
342     if mode == 'save':
343         plt.savefig(self.saveFile + '-XvectorFig.png')
344     elif mode == 'show':
345         plt.show()
346     else:
347         plt.savefig(self.saveFile + '-XVectorFig.png')
348         plt.show()
349
350     #Plot Y Vector
351     fig = plt.figure()
352     ax = fig.add_subplot(1, 1, 1)
353     ax.plot(time, self.y)
354
355     if title:
356         fig.suptitle(self.ID + ' -- ' + 'Y Vector')
357
358     if mode == 'save':
359         plt.savefig(self.saveFile + '-YVectorFig.png')
360     elif mode == 'show':
361         plt.show()
362     else:
363         plt.savefig(self.saveFile + '-YVectorFig.png')
364         plt.show()
365
366     #Plot Z Vector
367     fig = plt.figure()
368     ax = fig.add_subplot(1, 1, 1)
369     ax.plot(time, self.z)
370
371     if title:
372         fig.suptitle(self.ID + ' -- ' + 'Z Vector')
373

```

```

374         if mode == 'save':
375             plt.savefig(self.saveFile + '-ZvectorFig.png')
376         elif mode == 'show':
377             plt.show()
378         else:
379             plt.savefig(self.saveFile + '-ZvectorFig.png')
380             plt.show()
381
382     def createTimeAvgFig(self, xTimeAvg, yTimeAvg, zTimeAvg, time, mode='save',
legend=True, title=True):
383         fTime = self.formatTime(time)
384
385         fig = plt.figure()
386         ax = fig.add_subplot(1, 1, 1)
387         #plt.yscale('log')
388
389         if title:
390             fig.suptitle(self.ID + ' -- ' + 'Time Average')
391
392         ax.plot(fTime, xTimeAvg, label='X')
393         ax.plot(fTime, yTimeAvg, label='Y')
394         ax.plot(fTime, zTimeAvg, label='Z')
395
396         if legend:
397             ax.legend()
398
399         if mode == 'save':
400             plt.savefig(self.saveFile + '-timeAvgFig.png')
401         elif mode == 'show':
402             plt.show()
403         else:
404             plt.savefig(self.saveFile + '-timeAvgFig.png')
405             plt.show()
406
407     def createMagFig(self, magnitude, time, mode='save', title=True):
408         fTime = self.formatTime(time)
409
410         fig = plt.figure()
411         ax = fig.add_subplot(1, 1, 1)
412         plt.yscale('log')
413
414         if title:
415             fig.suptitle(self.ID + ' -- ' + 'Magnitude')
416
417         ax.plot(fTime, magnitude)
418
419         if mode == 'save':
420             plt.savefig(self.saveFile + '-timeMagFig.png')
421         elif mode == 'show':
422             plt.show()
423         else:
424             plt.savefig(self.saveFile + '-timeMagFig.png')
425             plt.show()
426

```

Supplementary Program 7 – Heatmap Generation (Python)

```
1 from simulation import Sim
2 from accelDataProcessing import AccelDataProcessor
3 from pathVisualization import PathVisualization
4 import numpy as np
5
6 class HeatMapGen:
7     innerMin = 0.125
8     innerMax = 4
9     outerMin = 0.125
10    outerMax = 4
11
12    vectorSim = Sim()
13    magList = []
14    distList = []
15    mapOut = open('C:\\MAP_PATH\\MAP_FILENAME', 'w+') #SET MAP PATH AND FILE NAME
16    rankOut = open('C:\\RANK_PATH\\RANK_FILENAME', 'w+') #SET RANK PATH AND FILE NAME
17    for inV in range(int(innerMin * 8), int(innerMax * 8) + 1):
18        magRow = []
19        distRow = []
20        for outV in range(int(outerMin * 8), int(outerMax * 8) + 1):
21            timeArr, xArr, yArr, zArr = vectorSim.gVectorData(0, 10000, float(inV)/8,
float(outV)/8)
22            process = AccelDataProcessor(xArr, yArr, zArr, timeArr)
23            xTimAvg, yTimAvg, zTimAvg = process.getTimeAvg()
24
25            magArr = []
26            for i in range(len(xTimAvg)):
27                curX = xTimAvg[i]
28                curY = yTimAvg[i]
29                curZ = zTimAvg[i]
30                curMag = (curX**2 + curY**2 + curZ**2)**0.5
31                magArr.append(curMag)
32
33            magRow.append(np.mean(magArr[5000:6000]))
34
35            id = str(inV) + ":" + str(outV)
36            pathVis = PathVisualization(id, xArr, yArr, zArr)
37            distRow.append(pathVis.getDistribution())
38            print(id)
39            magList.append(magRow)
40            distList.append(distRow)
41            magRow = []
42            distRow = []
43    print(magList)
44    print(distList)
45
46    mapOut.write("Magnitude\n")
47    for rowI in range(len(magList)):
48        for colI in range(len(magList)):
49            mapOut.write(str(magList[rowI][colI]) + "\t")
50        mapOut.write("\n")
51
52    mapOut.write("\nDistribution\n")
53    for rowI in range(len(distList)):
54        for colI in range(len(distList)):
55            mapOut.write(str(distList[rowI][colI]) + "\t")
56        mapOut.write("\n")
57
58    mapOut.close()
```



```

59
60 rankOut.write("Combination\tMagnitude\Distribution\n")
61 for rowI in range(len(magList)):
62     for colI in range(len(magList)):
63         inV = float(rowI + 1) / 8
64         outV = float(colI + 1) / 8
65         id = str(inV) + ":" + str(outV) + "\t"
66         rankOut.write(id + str(magList[rowI][colI]) + "\t" + str(distList[rowI]
[colI]) + "\n")
67     rankOut.close()
68

```

Supplementary Program 8 – Temperature Sensor (Python)

```
1 import os
2 import glob
3 import time
4 import numpy as np
5
6 openTime = time.strftime('%b/%d/%Y-%H:%M:%S')
7
8 os.system('modprobe w1-gpio')
9 os.system('modprobe w1-therm')
10
11 base_dir = '/sys/bus/w1/devices/'
12 device_folder = glob.glob(base_dir + '28*')[0]
13 device_file = device_folder + '/w1_slave'
14
15 curData = []
16
17 class TEMP:
18     def read_temp_raw():
19         f = open(device_file, 'r')
20         lines = f.readlines()
21         f.close()
22         return lines
23
24     def read_temp():
25         lines = read_temp_raw()
26         while lines[0].strip()[-3:] != 'YES':
27             time.sleep(0.2)
28             lines = read_temp_raw()
29         equals_pos = lines[1].find('t=')
30         if equals_pos != -1:
31             temp_string = lines[1][equals_pos+2:]
32             temp_c = float(temp_string) / 1000.0
33             temp_f = temp_c * 9.0 / 5.0 + 32.0
34             return temp_c
35
36     def startSession():
37         d = open('incTemp1.txt', 'a+');
38         d.write('\n' + 'Session Start: ' + str(openTime) + '\n\n')
39
40     def curTime():
41         return time.strftime('%b/%d/%Y-%H:%M:%S')
42
43     def getTemp():
44         if len(curData) < 1:
45             startSession()
46         curData.append(read_temp())
47         d.write(str(read_temp()) + '\t' + curTime() + '\n')
```